



Application Data Sheet




High Performance Liquid Chromatography

No. 39

High-Speed Analysis of Geniposide

Here we conducted high-speed analysis of geniposide, one of the iridoid glycosides contained in the natural medicine gardenia (gardenia fruit).



Analysis of Geniposide Found in Traditional Chinese Medicine

[Sample Preparation]

To produce a mixture of 10 g/L, finely powdered tablet was mixed with 50% methanol solution, and then filtered through a 0.45 µm membrane filter.

Analytical Conditions

Instrument : Prominence UFLC system

Column : Shim-pack XR-ODS (50 mm × 3.0 mm i.d.)

Mobile Phase : Water / Acetonitrile / phosphoric acid = 900 / 100 / 1 (v/v)

Flow Rate : 1.2 mL/min Column Temperature : 35 $^{\circ}$ C

Detection : Absorbance at 240 nm with Semi-micro flow cell

Sample Volume : 4 µL

